B.Sc. (Part—III) Semester—VI Examination BIOTECHNOLOGY (R/V)

(Plant Biotechnology)

	(Fiant Die	necumology)	
Time: Three	Hours]		[Maximum Marks: 80
Note :—(1) (2)	ALL questions are compulsory. Draw neat and labelled diagram	s wherever necessary.	
1. (A) Fill	in the blanks :		2
(i)	Specialised pores on leaf meant	for transpiration are ca	lled as
(ii)	is the fruit ripening hor	mone.	
(iii)	The products of Mutation are ca	alled as	
(iv)	Haploids are produced by	culture.	
(B) Mul	tiple choice questions :		. 2
(i)	Indole-3 acetic acid is an examp	ole of:	
	(a) Auxin	(b) Cytokinin	
	(c) Abscisic Acid	(d) Gibberalic Ad	id
(ii)	Endosperm culture of an Angios	spermic plant produce :	
	(a) Haploids	(b) Diploids	
	(c) Triploids	(d) Tetraploid	
(iii)	Enzymes required for protoplast	isolation:	
	(a) Lypase	(b) Isomerase	
	(c) Cellulase and pectinase	(d) Proteases	
(iv)	Which particles are not used in		
	(a) Gold	(b) Silver	
	(c) Tungsten	(d) Carbon	
VTM—13441		1	(Contd.)

www.sgbauonline.com

	(C)	Answer in ONE sentence each :	4
		(i) Define Totipotency.	
		(ii) What is Protoplast ?	
		(iii) What is meant by organogenesis?	-
		(iv) What is phototropism?	
2.	Des	cribe :	
	(a)	Apical dominance	4
	(b)	Growth curve	4
	(c)	Geotropism	4
		OR.	
	(d)	Long day plants	4
	(e)	Structure of stomata	4
	(f)	Phototropism.	4
3.	Exp	lain in detail the physiological effects of cytokinins.	12
		OR	
	Exp	lain the mechanism of action of Auxin.	12
4.	Exp	lain in brief:—	
	(a)	Composition of M.S. Medium.	4
	(b)	Practical applications of Tissue Culture.	4
	(c)	Management of Tissue Culture Laboratory.	4
		OR	
	(d)	Role of Laminar Air Flow in Tissue Culture.	4
	(e)	Practical applications of Organ Culture.	4
	(f)	Sterilization techniques in Tissue Culture.	4
V Th	i 13	141	(0
* I !V	1,3	141 2	(Contd.)

www.sgbauonline.com

5.	Exp	olain :—	
	(a)	Ovule culture	4
	(b)	Applications of Somaclonal variation	4
	(c)	Micropropagation.	4
		OR	
	(d)	Embryo rescue	4
	(e)	Clonal multiplication	4
	(f)	Anther culture.	4
6.	Des	cribe :—	
	(a)	Passive transport across cell membrane.	4
	(b)	Protoplast isolation (enzymatic method),	4
	(c)	Single cell suspension culture.	4
		OR	
	(d)	Active transport across cell membrane.	4
	(e)	Role of single cell suspension culture in selection of variants.	4
	(f)	Protoplast regeneration.	4
7.	Des	cribe Agrobacterium mediated gene transfer in detail.	12
		OR	
	Exp	lain somatic hybridization and its applications in detail.	12

www.sgbauonline.com